



Oxford Instruments Superconductivity

UK
Tubney Woods, Abingdon
Oxfordshire OX13 5QX
Tel: +44 (0)1865 393 200
Fax: +44 (0)1865 393 333

China
Room 8418a, Building A
Jiahua Business Center
808 Hongqiao Road
Shanghai 200030
Tel: +86 21 64867993
Fax: +86 21 64861781

Germany
Otto-von Guericke Ring 10
D-65205 Wiesbaden
Tel: +49 6122 937 171
Fax: +49 6122 937 175

Italy
Via Leone Tolstoj 86
20098 San Giuliano
Milanese
Milan
Tel: +39 02 98 2531
Fax: +39 02 98 24407

Japan
Haseman Building
2-11-6 Tomioka, Koto-ku
Tokyo 135-0047
Tel: +81 03 5245 3261
Fax: +81 03 5245 4472

Spain
Avda Mata Piñonera, 2
28700 San Sebastian de los Reyes
Madrid
Tel: +34 91 659 0740
Fax: +34 91 654 6794

U.S.A.
130A Baker Ave. Ext.
Concord, MA 01742-2121
Tel: +1 978 369 9933
Fax: +1 978 369 6616

E-mail:
superconductivity@oxinst.co.uk

Visit our web site at
www.oxford-instruments.com

OXFORD
INSTRUMENTS

This publication is the copyright of Oxford Instruments Superconductivity Limited and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations

© Oxford Instruments Superconductivity Limited, 2003. All rights reserved.

Reference No: HCC OI627 02/03

UK sites of Oxford Instruments Superconductivity Limited, operate a Quality Management System approved to the requirements of BS EN ISO 9001.



Certificate No Q4118

Sample in vacuum 1 K pot-free ³He refrigerator

The HelioxVT is specifically designed to fit into Oxford Instruments Superconductivity's variable temperature inserts (VTI) with 30 mm access. Operated in the cold gas environment of a variable temperature insert, this removes the need for the 1 K pot on the insert. This eliminates the presence of liquid helium in the sample horizontal plane. These unique features offer two key benefits:

- Extension of sample temperature range to below 280 mK of existing or new VTIs
- Ideal for neutron or X-ray scattering experiments




Figure 1: The HelioxVT system

Components

A complete system consists of:

- An HelioxVT insert (Figure 1)
- An ITC⁵⁰³ temperature controller to operate the ³He refrigerator
- A variable temperature insert (VTI) with pump and ITC⁵⁰³

Features and Benefits

- **Used in a VTI environment** – extends temperature range to temperatures below 280 mK.
- **No liquid helium in beam plane** – absence of cryogen in the sample horizontal plane enables scattering experiments (neutron or X-ray) to be undertaken without disruptive beam dissipation or reflection from Helium
- **Small diameter insert** – specially designed 30 mm insert permits a very compact and small footprint system.
- **Wide temperature range** – temperature can be controlled from 280 mK up to 6 K (or higher with high temperature option)
- **Seamless integration** – with Oxford Instruments Superconductivity designed VTI and split pair type superconducting magnets.



OXFORD
INSTRUMENTS

Operation

Figure 2 represents a schematic of the HelioxVT running in a VTI. The HelioxVT insert can be treated like any sample rod for a variable temperature insert. Once loaded, the insert is cooled by the VTI from 300 K down to around 10 K using exchange gas. The ^3He gas contained in a small dump sitting on top of the insert is then condensed at around 1.5 K. Once the ^3He pot has reached a stable temperature and condensation is completed, the adsorption pump will start to cool the ^3He pot and experimental set-up to below 300 mK. The condensation and the cool down time typically require less than 1 hour.

Performance

Despite its exceptional compactness, the HelioxVT can maintain temperatures of below 280 mK for more than 40 hours (with no applied heat load). In addition, the integrated sliding seal enables complete ease of use and flexibility of operation. Available cooling power of more than $50 \mu\text{W}$ at 350 mK can be maintained for over 4 hours, whilst typical regeneration time is around 40 minutes.

Experimental access

The insert is equipped with various diagnostic sensors, including a carbon sensor on the sorption pump and uncalibrated RuO_2 sensor on the 1.5 K condensing stage and ^3He pot. A universal (50 mm) connector is fitted to the top of the insert. This connector is equipped as standard with two 10-way Fischer connector (diagnostic and customer wiring) 5 twisted pairs, 0.1 mm diameter constantan wires terminated at the ^3He pot.

HelioxVT

Standard Specification		Description
HELV30SYS	Base temperature	<280 mK for >24 hrs (no applied heat load)
	Cooling power	<350 mK for >4 hours with $50 \mu\text{W}$ applied heat load
	Temperature range	Base temperature to 6 K
	Temperature stability	Below 1.2 K: ± 3 mK
	Thermometers	Uncalibrated RuO_2 sensors on 1 K stage and ^3He pot. Carbon resistor on sorption pump.
	^3He regeneration time	40 min
	Hold time at base temperature	>24 hours
	Sample access	Via IVC; Inner diameter 25 mm
	Experimental wiring	10-way loom wired in twisted constantan pairs to 10-pin socket at ^3He pot

High Temperature Upgrade		HelioxVT extended operating temperature range
HELHT	Thermometers	RhFe or Cernox™ sensor to ^3He pot or probe for operation above 6 K.
	Temperature range	Up to 80 K
	Temperature stability	± 100 mK at $T > 1.5$ K

Standard Options	
SHM	Sample holder for high field magnet
SIVCT	Spare IVC tail, to suit magnet if required




Figure 2: Operating principle of HelioxVT