

```
In [1]: # this guide requires python 3.6
import sys
sys.version
```

```
Out[1]: '3.6.7 | packaged by conda-forge | (default, Nov 6 2019, 16:19:42) \n[GCC
7.3.0]'
```

```
In [2]: import numpy as np
import h5py
```

```
In [3]: # open the file read-only, with error handling
try:
    h5 = h5py.File("pshell-20200824-191714-ManipulatorScan.h5", "r")
except OSError:
    print("error opening file")
```

```
In [4]: # just checkin'
print(h5.file.filename)
```

```
pshell-20200824-191714-ManipulatorScan.h5
```

```
In [5]: # what's inside the file at top level?
h5.keys()
```

```
Out[5]: <KeysViewHDF5 ['__DATA_TYPES__', 'general', 'logs', 'scan 1']>
```

```
In [6]: # the general group contains the experiment metadata
general = h5['general']
general.keys()
```

```
Out[6]: <KeysViewHDF5 ['authors', 'pgroup', 'proposal', 'proposer', 'sample']>
```

```
In [7]: # authors is an array, the other ones are singletons
# call the decode method to convert to string
print([v.decode() for v in general['authors']])
print(general['pgroup'][()].decode())
print(general['sample'][()].decode())
```

```
['Matthias Muntwiler']
p17274
Cu(111)
```


```
In [8]: # most files contain just one scan group.
# XPS spectra, e.g., can contain multiple scans.
#
# what's inside a scan group?
h5['scan 1'].keys()
```

```
Out[8]: <KeysViewHDF5 ['Counts', 'ImageAngleDistribution', 'ImageEnergyDistribution',
'MachineCurrent', 'ManipulatorPhi', 'RefCurrent', 'SampleCurrent', 'Scienta image', 'Scienta spectrum', 'attrs', 'meta']>
```

```
In [9]: import matplotlib.pyplot as plt
```

```
In [10]: # the datasets (all items except 'attrs' and 'meta') contain the measured data
# they are represented as numpy datasets and can be plotted directly
plt.plot(h5['scan 1']['ManipulatorPhi'], h5['scan 1']['Counts'])
```

```
Out[10]: [<matplotlib.lines.Line2D at 0x7ff16cb72a90>]
```



```
In [11]: # the attrs group contains beamline status
h5['scan 1']['attrs'].keys()
```

```
Out[11]: <KeysViewHDF5 ['AcquisitionMode', 'AnalyserSlit', 'BeamlinePressure', 'ChamberPressure', 'DetectorMode', 'DynamicBumpYAngle', 'DynamicBumpYOffset', 'ElementSet', 'EnergyMode', 'ExcitationEnergy', 'ExitSlit', 'FocusWaterTemp', 'FocusXRot', 'FocusYRot', 'FocusYTrans', 'FocusZRot', 'FocusZTrans', 'FrontendHCenter', 'FrontendHSize', 'FrontendVCenter', 'FrontendVSize', 'LensMode', 'MachineBumpXAngle', 'MachineBumpXOffset', 'MachineBumpYAngle', 'MachineBumpYOffset', 'MachineCurrent', 'ManipulatorCoolFlow', 'ManipulatorCoolFlowSet', 'ManipulatorPhi', 'ManipulatorTempA', 'ManipulatorTempB', 'ManipulatorTheta', 'ManipulatorTilt', 'ManipulatorX', 'ManipulatorY', 'ManipulatorZ', 'MonoApertureMode', 'MonoBeta', 'MonoBladeDown', 'MonoBladeUp', 'MonoCff', 'MonoEnergy', 'MonoGrating', 'MonoHCenter', 'MonoHSize', 'MonoTheta', 'MonoVCenter', 'MonoVSize', 'NumIterations', 'PassEnergy', 'RefCurrent', 'RefCurrentGain', 'RefocusHCenter', 'RefocusHSize', 'RefocusVCenter', 'RefocusVSize', 'RefocusXRot', 'RefocusYRot', 'RefocusYTrans', 'RefocusZRot', 'RefocusZTrans', 'SampleCurrent', 'SampleCurrentGain', 'ScientaCenterEnergy', 'ScientaChannelBegin', 'ScientaChannelEnd', 'ScientaDwellTime', 'ScientaHighEnergy', 'ScientaLowEnergy', 'ScientaNumChannels', 'ScientaNumSlices', 'ScientaSliceBegin', 'ScientaSliceEnd', 'StepSize']>
```

```
In [12]: # create a shortcut to the first scan group
# scan1 and h5['scan 1'] can be used interchangeably
scan1 = h5['scan 1']
print(scan1)
```

```
<HDF5 group "/scan 1" (11 members)>
```

```
In [13]: # a scan group also contains direct 'attributes'
# note that 'attributes' is used with several different meanings - unfortunately:
# - the scan group contains another group called 'attrs'
# - the scan group has 'attributes', which are represented in h5py by the 'attrs' attribute :-o
#
scan1.attrs.keys()
```

```
Out[13]: <KeysViewHDF5 ['Dimensions', 'End', 'Passes', 'PlotTypes', 'Readables', 'Start', 'Steps', 'Writables']>
```

```
In [14]: # number of scan dimensions
print(scan1.attrs['Dimensions'])
# number of scan steps
print(scan1.attrs['Steps'])
# number of iterations (XPS spectra only)
try:
    print(scan1.attrs['Iterations'])
except KeyError:
    print("(no iterations)")
# number of passes
print(scan1.attrs['Passes'])
# start and end time in microseconds
print((scan1.attrs['End'] - scan1.attrs['Start']) / 1000 / 60)

1
[10]
(no iterations)
1
7.917883333333333
```

```
In [15]: # list of scan positioners
scan1.attrs['Writables']
```

```
Out[15]: array([b'ManipulatorPhi'], dtype='<S14')>
```

```
In [16]: # remember to cast to string
scan1.attrs['Writables'][0].decode()
```

```
Out[16]: 'ManipulatorPhi'
```

```
In [17]: # list of detectors
scan1.attrs['Readables']
```

```
Out[17]: array([b'Counts', b'MachineCurrent', b'SampleCurrent', b'RefCurrent',
                b'Scienta spectrum', b'ImageEnergyDistribution',
                b'ImageAngleDistribution', b'Scienta image'], dtype='<S23')>
```

```
In [18]: # note - these are the names of the datasets in the scan group!
scan1.keys()
```

```
Out[18]: <KeysViewHDF5 ['Counts', 'ImageAngleDistribution', 'ImageEnergyDistribution',
                    'MachineCurrent', 'ManipulatorPhi', 'RefCurrent', 'SampleCurrent', 'Scienta image',
                    'Scienta spectrum', 'attrs', 'meta']>
```

```
In [ ]:
```