

```

sics 'user Muster
proposolid 2020001';

use_new_pswEEP;

sics 'radcol rc1';      # specify the current radial collimator installed
#sics 'radial start';  # starts the radial collimator
sics 'radial check 1'; # make sure sics checks the status of
                        # radial collimator,
                        # and restarts it, if it is stopped.

sample_changer 'lt5';  # we choose 'lt5', other options are 'rt', 'lt'
sing_current 1000;     # approximate SINQ current to have good timing
time_per_sweep '1 hour';
                        # time per one datanumber/file. We call it sweep.

sample_rotation 'on';
beam_reduction_horizontal 12;
beam_reduction_vertical 50;
                        # for small height of the sample H the vertical
                        # can be ~ H+10. All in [mm]

lambda 1.5;           # neutron wavelength in Angstroms
resolution 'HI';     # other choices are MR, HR
# lambda 1.9;
# set_sweep_range '3.55,0.15,4';
time_per_sweep '0.5 hour';

ttol 15;              # temperature tolerance

# hcount 100, 'empty_IRI4';
# lambda 2.45;
# hcount 1, 'Li-biPy_n1', 4, 300;
# hcount 10, 'NaKCu30(SO4)3_V8x25', 3, 300;
# hcount 19, 'Li-biPy_n1', 4, 300;

# hcount -10, 'cool_Albert_V8x25', 3, 1;
# hcount 10, 'NaKCu30(SO4)3_V8x25', 3, 1;

ttol 2;
lambda 1.9;
hcount 3, 'Yb0p667Mn0p333Mn03_V8x10', 1, 1;
                        # main count command
                        # explained below
# hcount 3,          'Yb0p667Mn0p333Mn03_V8x10',          1,          1;
#           |              |              |              |
#           v              v              v              v
#           3 sweeps      full sample name      position 1      temperature 1 in [K]
#                               in sample changer

```

```

for $t (10, 20, 30, 40, 50, 60, 75, 90, 130) {
    # the list of temperatures to be
    # measured

    ttol 0.5 if $t>1;
    ttol 1 if $t>30;
    ttol 2 if $t>50;
    ttol 5 if $t>75;

    # the above 4 lines set the tolerance
    # depending on the set point, to avoid
    # long waiting times - temperature of the sample
    # reaches the set point exponentially slow...

    hcount 4, 'Yb0p667Mn0p333MnO3_V8x10',1,$t; # in position 1
# 6K, 15K, 25K, 35K, 45K, 55K
}
ttol 10;
for $t (qw/130 170 210/) {
    hcount 4, 'Tm0p7Mn0p3MnO3_V8x15',2,$t;
}
ttol 2;
hcount -10, 'cool_Tb0.6Pr0.4MnO3',5,1;
lambda 2.45;
hcount 99, 'Li-biPy_n1',4,1;

```