


## 1.) Load sample into Prepchamber


*How to load a sample using the ELMITEC load lock for transfer into the PEEM.*

### Location

- The load lock is located below the prepchamber on the left-hand side of the PEEM.



- The forevacuum buffer is in the rack WSLA 13.0.55



### 1.1.) Venting the load-lock

*Venting the load-lock to insert the new sample.*

- Turn N2 gas ON (behind hutch) and check for flow.
- Set Buffer controller to P2 = ACTIVE.
- Check that valve separating load-lock / prepchamber and prepchamber / mainchamber are CLOSED.
- Turn TMP controller (DCU) Power = OFF
- Wait until load-lock is vented.  
Vent valve opens once speed is below 1350/min, after ~3 min.

- Open CF40 flange
- Turn N2 gas OFF (behind hatch).

**1.2.)Install new sample**

*Install new sample and prepare for pump-down.*

- Close CF40 flange (two M6 screws and VITON gasket).
- Set Buffer controller to P1 + P2 = ACTIVE.  
Wait 50sec until valve opens (Vent 1 lights up).
- Turn TMP controller (DCU) Power = ON  
Check that pressure goes down ( $10^{-2}$  mbar after ~ 1min).
- When  $P_{LOAD} < 2 \cdot 10^{-6}$  mbar transfer sample.

**Known problems:**

- *Sometimes the vent valve of the TMP does not close. Turn TMP controller OFF / ON*