



- To prevent cold-welding between the piston and the pressing pad, put a piece of paper with lubricant between them
- Combine the bottom bolt with the pressing pad.
- Apply some molycote grease on the threads of the bottom bolt
- Close the bottom bolt.

**Now flip the cell and assemble the top part.**

- Place pressure gauge (Indium disc)
- Place sample (Ideally, a few cylinders of combined 10-12 mm height and the required diameter)
- Fill the remaining of the volume with Daphne 7373 (or Daphne 7474) oil.
- Remove some Daphne oil (about xx mm, using the pre-gauged syringe)
- Place the second mushroom.
- Put the sealing ring + disc on top of the mushroom. (sharp edge of the ring towards the center)
- Add one drop of lubricant (molybdenum disulfite +daphne oil mixture)
- Put the piston (the longer one) in
- To prevent cold-welding add a piece of paper with lubricant
- Combine the top bolt with the pressing pad.
- Apply some molycote on the threads of the top bolt
- Close the top bolt

**The cell is now assembled but not sealed yet.**

- To seal the cell, fix the cell body to the table and tighten the bottom bolt.
- Once the bottom bolt is tight, then switch to the top bolt and tighten it by hand.
- The sealing procedure, depending on how strong one closes it already gives about 5 kbar at room temperature, so a few kbar at low temperatures.

**The cell is now sealed and ready for high pressure application.**

### 3. Application of pressure

- Place the cell into a cylinder (#1) as shown in the figure on the right:
- The cylinder already centers the cell.
- Put a pressure application piston (#2) into the hole inside the top bolt

**Now everything is assembled, the press will soon be turned on**

- Check the conversion factor to obtain the right value of pressure in the cell from the manometer.
- Close all the potentially dangerous solid angle by plexi glass.
- Move security valve to 'retract' position
- Turn on the press. The press will engage with the pressure application piston and apply a small starting pressure of the order of 0.5 kbar.
- Move the security valve to 'hold' position
- Fix the distance measurement device (#3).
- Start the recording
- Make sure that the needle valve is closed.
- Move the security valve to 'advance' position
- Slowly open the needle valve. Monitor both the pressure on the screen and on the manometer. Aim at a rate of 0.2 bar per second on the manometer for a slow and continuous pressure application.
- Once the desired pressure is reached, move the security valve to 'hold position'
- Now use a wrench to transfer the force from the press to the top bolt. Be careful and wear gloves.
- Look at the manometer and try to reduce the reading by about 5 bar.
- Once it is tightened, switch the security valve to 'retract' position and allow the pressure to be slowly removed from the press.
- After the pressure has dropped to 1 kbar, the press can be switched off. Then open the needle valve fully for the head to move away from the cell.

**Now the cell is pressurized. Leave it for 10 minutes for moderate pressures and 30 minutes for pressures larger than 75% of the maximum pressure.**



#### 4. Further manipulation

Depending on what needs to be done next (load to the spectrometer or to the ACS measurement), different steps will be taken. But in all occasions, try to minimize the exposure of the cell to outside. Namely – use specialized ammo-boxes to transport the cell and perform all the handling behind the plexiglass shields.